


WHY Lineassay?



- easy processing
- time and cost saving multi-parameter assay
- high specificity due to highly purified natural antigens
- good diagnostic performance by simultaneous determination of several autoantibodies


TEST PRINCIPLE AND EVALUATION

- highly purified lipid antigens are immobilized on the test strip
- specific autoantibodies of patients sera bind to the immobilized antigens
- bound autoantibodies are detected by the secondary antibody and result in an enzymatic colour reaction during substrate incubation



- for evaluation an interpretation template is used
- the template shows the cut-off-value for all parameters (see below)

attomol® Anti-Phospholipid IgG/M Lineassay 3 (REF1176)



attomol® Anti-Ganglioside IgG/M Lineassay 6 (REF1189)

